

INSTRUCTION MANUAL



Equipment Type: 300 mm (~12-inch) Manual Abrasive Cut-off

Machine

Model: **MEGA-T300**

Electrical Requirements: 208Y / 380 Volts (3-phase)

Optional 480V 3-phase)

Frequency: 50/60 Hz

Motor Horsepower: 5.5 hp (4 KW) (3-phase)

Manual Revision Date: April 21, 2014 This instruction manual is provided with each piece of

delivered equipment.

Please read this instruction manual carefully and follow all installation, operating and safety guidelines.



INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

WARRANTY

Terms and Conditions applying to all PACE Technologies Products

1. LIMITED WARRANTY AND DISCLAIMER:

PACE Technologies Products are warranted for one year from the purchase date to be free from defects in material and workmanship under correct use, normal operating conditions, and proper application. PACE Technologies obligation under this warranty shall be limited to the repair or exchange, at PACE Technologies option, of any PACE Technologies Product or part which proves to be defective as provided herein. PACE Technologies reserves the right to either inspect the product at Buyer's location or require it to be returned to the factory for inspection. Buyer is responsible for freight to and from factory on all warranty claims. The above warranty does not extend to goods damaged or subjected to accident, abuse or misuse after release from PACE Technologies warehouse, nor goods altered or repaired by anyone other than specifically authorized PACE Technologies representatives. PACE Technologies shall not in any way be responsible for the consequences of any alteration, modification or misuse unless previously approved in writing by an officer of PACE Technologies.

PACE TECHNOLOGIES MAKES NO EXPRESS WARRANTIES OTHER THAN THOSE WHICH ARE SPECIFICALLY DESCRIBED HEREIN. Any description of the goods sold hereunder, including any reference to Buyer's specifications and any description in catalogs, circulars and other written material published by PACE Technologies, is the sole purpose of identifying such goods and shall not create an express warranty that the goods shall conform to such description.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED. THERE ARE NO IMPLIED WARRANTIES OF MECHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. THIS WARRANTY STATES PACE TECHNOLOGIES ENTIRE AND EXCLUSIVE LIABILITY AND BUYER'S EXCLUSIVE REMEDY FOR ANY CLAIM FOR DAMAGES IN CONNECTIONS WITH PACE TECHNOLOGIES PRODUCTS. PACE TECHNOLOGIES WILL IN NO EVENT BE LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES WHATSOEVER, NOR FOR ANY SUM IN EXCESS OF THE PURCHASE PRICE.

2. LIABILITY CAP:

PACE Technologies maximum aggregate liability for loss and damage arising under, resulting from or in connection with the supply or use of the Equipment and Consumables provided under this purchase, or from the performance or breach of any obligation (s) imposed hereunder, whether such liability arises from any one or more claims or actions for breach of contract, tort, (including negligence), delayed completion, warranty, indemnity, strict liability or otherwise, unless otherwise limited by the terms hereof, shall be limited to one hundred percent (100%) of the purchase price.

3. DELIVERY:

Customer assumes and shall bear the risk of all loss or damage to the Products from every cause whatsoever, whether or not insured, and title to such Products shall pass to Customer upon PACE Technologies delivery of the Products to the common carrier of Pace Technologies choice, or the carrier specified in writing by Customer, for shipment to Customer. Any claims for breakage, loss, delay, or damage shall be made to the carrier by the Customer and Pace Technologies will render customer reasonable assistance in prosecuting such claims.



INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

4. ACCEPTANCE:

Customer shall inspect the Products promptly upon receipt of delivery. Unless customer objects in writing within five (5) business days thereafter, customer shall be deemed to have accepted the Products. All claims for damages, errors, or shortage in Products delivered shall be made by Customer in writing within such five (5) business days period. Failure to make any claim timely shall constitute acceptance of the Products.

Any returns for credit may be subject to 25% restocking fee and must be returned within ten (10) business days from initial delivery.

5. PAYMENT:

Customer agrees to provide timely payment for the Products in accordance with the terms of payment set forth on the reverse side hereof or in any proposal submitted herewith. If any payment is not paid on or before its due date, Customer shall pay interest on such late payment from the due date until paid at the lesser of 12% per annum or the maximum rate allowed by law.

6. MISCELLANEOUS PROVISIONS:

This agreement has been made in and shall be governed by the laws of the State of Arizona. These terms and conditions and the description of the Products on the reverse side hereof or in any proposal submitted herewith constitute the entire agreement and understanding of the parties with respect to this sale and supersede all prior and contemporaneous agreements or understandings, inducements or representations, expressed or implied, written or oral, between the parties with respect hereto. Any term or provision of this Agreement may be amended, and any observance of any term of this Agreement may be waived, only by a writing signed by the party to be bounds. The waiver by a party of any breach shall not be deemed to constitute a waiver of any other breach. Should suit be brought on this Agreement, the prevailing party shall be entitled to recover its reasonable attorneys' fees and other costs of suit including costs and attorneys' fees incurred on appeal or in collection of any judgment.

INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

Contents

Nan	neplate	PAGE i
Waı	ranty	ii
1.0	Product Description	1
2.0	Shipping, Unpacking and Installation	6
3.0	Safety Guidelines	9
4.0	Start-Up and Operation	16
5.0	Maintenance	17
6.0	Trouble Shooting	18
7.0	List of Spare Parts	19
8.0	Electrical Drawings	41

INSTRUCTION MANUAL


3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

1.0 Product Description


1.1 General Description

The **MEGA-T300** is a manual wheel feed or table feed abrasive cut-off machine for cutting materials ranging from soft aluminum metals to hardened tool steels. It is ideal for the metallographic laboratory, as well as for small industrial or production applications.

The **MEGA-T300** is very robust and durable with its cast aluminum alloy and stainless steel construction. Featuring a corrosion-free T-slot table, the **MEGA-T300** is a very versatile table feed bench top model metallographic cutter.



INSTRUCTION MANUAL



INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

1.2 Technical Specifications

Electrical specifications: 208Y (phase-phase or 120/120/120 to

ground), 380V, 400, 480V 3-phase

optional)

Motor power: 4.5 hp (3 kW) - 3-phase to ground

Cut-off wheels: 300 mm (~12-inch) diameter

Wheel arbor: 32 mm (~1.25 inch) diameter

Maximum sample diameter: 4-inch (100 mm)

Weight: Approx. 350 lbs (160 kg)

Dimensions (WxHxD): Approx. 36" x 34" x 26" (hood closed)

(90 cm x 86 cm x 65 cm)

Approx. 36" x 34" x 32" (hood open)

(90 cm x 86 cm x 82 cm)

Table dimensions (WxD): 12" x 9"

(300 mm x 230 mm)

Cutting chamber: Cast aluminum block construction

Hood and body Fabricated steel with shock resistant

plastic window

Working temperature: 32° - 100°F (0 - 40°C)

Shipping temperature: 32° - 100°F (0 - 40°C)

Storage temperature: 32° - 100°F (0 - 40°C)

Recirculation system (included): 15 gallons (60 liters)

 $C \in$

EU Directives: Machinery directive 2006/42/EC

Low Voltage Directive 2006/95/EC

Electromagnetic Compatibility directive 2004/108/EC


EU Harmonized Standards: EN ISO 1200:2010


EN 61010-1:2010 EN 61326-1:2006


INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

1.3 Mechanical Schematic







Note: Installation of the *MEGA-T300* should be on a flat sturdy surface, with easy access to electrical connections.

INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

1.4 Features



Powerful motor


The **MEGA-T300** is equipped with a powerful belt drive motor. In addition, the motor has an inductive brake for faster stopping of the blade.

Manual feed and table feed cutting

The **MEGA-T300** is an easy-to-use manual cutter with simple controls. Cutting can be accomplished with either a wheel feed or table feed.

The **MEGA-T300** also has a removable side port panel for sectioning longer samples.





Safety lock-out for protecting



INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

2.0 Unpacking, Shipping and Installation

2.1 Unpacking

Unit is delivered in a crate. Unpack and check for completeness of parts.

Measures WxHxD: Approx. 36" x 36" x 36"

Weight: Approximately 375 lbs.

2.2 Shipping

When moving crate, lift from bottom. To move MEGA-T300 use support bars provided with the machine.







The *MEGA-T300* is constructed of sensitive electronic and mechanical components. **Do not drop.**

Caution: Heavy equipment. Take care to avoid bodily injury.



INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

2.3 Installation

Install unit carefully! Improper installation voids warranty.

The **MEGA-T300** should be placed on a flat stable surface.

TO OPEN HOOD, CUTTER IS SHIPPED WITH THE LOCKOUT SWITCH TURNED TO UNLOCK. AFTER CONNECTING POWER, TURN SWITCH TO LOCK TO OPEN HOOD AND TO USE CUTTER

Connect coolant tank supply, drain and electrical connections.





CAUTION

Before installation or when there is no power to the machine, the hood cannot be opened with the interlocking safety device in the **LOCK** position. To open the hood without power hook-up, turn the interlocking device with key to the **UNLOCK** position (used for service and shipping).

Shipped in the Unlock Turn to Lock position to operate



Electrical connections—The standard saw is designed to operate at 208-220V 3-phase (120V-120V-120V-ground) (optional 380/400/480V)

Verify voltage on nameplate





Main Power switch


Verify the direction of rotation of the cut-off wheel. The wheel should t_____ as viewed from the front of the machine.

(Installation continued on next page)

INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

Recirculating tank connections



External coolant supply: Attach 3/4-inch tube between pump and cutter.

Drain: 1.5-inch (38 mm) tube.

Electrical connections: Connect six-foot electrical power cable to source.

Note: Inspect the operating voltage on the name plate.

Electrical connection for

external coolant supply: Power for recirculation system comes from the **MEGA-T300**,

(Installation continued on next page)



INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

3.0 Safety Guidelines

3.1 Warning Sign

This sign points to special safety features on the machine.

3.2 Safety Precautions

Careful attention to this instruction manual and the recommended safety guidelines is essential for the safe operation of the *MEGA-T300*.

Proper operator training is required for the safe operation of the *MEGA-T300*. Any unauthorized mechanical and electrical change, as well as improper operation, voids all warranty claims. All service issues need to be reported to the manufacturer / supplier.

- Before operating, cutting chamber hood must be closed. After cutting, the safety latch will not open for approximately 5 seconds after pressing the stop (red) button.
- Use only certified cut-off wheels from a professional supplier. Improper blades selection voids warranty. (For appropriate blade selection, refer to the Abrasive Blade Selection Guidelines Chart in Section 4.4)
- Disconnect power before opening the main unit.
- Replacement parts should be installed only by qualified personnel.
- Securely clamp the part / sample to the working table. During cutting, consider that the part may pinch and cause jamming of the cut-off wheel. Use the appropriate clamping devices to avoid this occurrence.
- Never start a cut under load.
- Make sure that the cut-off blade is rotating down and into the sample.

3.3 Emergency Statement

The MEGA-T300 abrasive cutter has been designed for cutting metallographic specimens. DO NOT CUT oversize samples. Always follow proper operational guidelines and avoid contact with moving parts, lubricants and abrasives. Seek appropriate medical care for cutting injuries.



INSTRUCTION MANUAL


3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

Recirculation tank

The **MEGA-T300** is equipped with a mobile three-chamber counter-current coolant tank. The first chamber contains a coarse strainer for removing larger debris.

Replacement of the coolant is accomplished by rolling out the tank and removing the water pump.

Note: Dispose of the old coolant in accordance with federal, state and local regulations.










Three Chamber tank with counter-current flow



Coarse strainer

Rinse line

After cutting, it is recommended that the cutting chamber be rinsed to prevent debris and corrosion build-up.




Water flow regulating valve



INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

3.4 Safety Tests



Examine and verify that the **MEGA-T300** safety devices and operating performance are in good working condition prior to use. The following safety checks are considered important:

Emergency stop switch

Test: Activate main switch and close hood;

depress emergency stop switch.

Proper

Response: Machine powers down.

Malfunction: Machine does not lose power.

Corrective If system does not power down,

measure: disconnect power supply cord and call

service technician.

Magnetic safety lock out switch

Test: Activate main switch and close hood;

turn cut-off wheel ON then OFF;

Try to open the hood.

Proper Hood does not open for approximately 5

Response: seconds.

Malfunction: Hood opens immediately.

Corrective If cut-off wheel does not power down,

measure: disconnect power supply cord and call

service technician.



Magnetic Safety switch

INSTRUCTION MANUAL


3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

Start-up and Operation

General

The **MEGA-T300** is a manual wheel and manual table feed cutter.

4.2 Control Panel



Cutting wheel ON/: Starts the cutting wheel.

Cutting wheel OFF/: Stops the cutting wheel (inductive brake engages to slow

> wheel). Hood cannot be opened until safety lock releases in approximately 5 seconds after pressing the stop button.

Coolant ON/OFF: Operates the coolant pump in Auto or Manual mode.

Light ON/OFF: Turns inside chamber light On or Off.

INSTRUCTION MANUAL

4.3 Wheel and Table Feed

4.3.1 Wheel feed cutting



4.3.2 Table feed cutting




Open wheel spinner

INSTRUCTION MANUAL

4.4 Changing abrasive cut-off wheels

- 1. Remove blade locking bolt (reverse threaded).
- 2. Position new 300 mm abrasive cut-off wheel into position.
 - Use only certified abrasive cut-off wheels.
- 3. Gently tighten blade locking bolt (Note: locking bolt is reverse threaded).





INSTRUCTION MANUAL

----- 🛦

ABRASIVE Cutting CONSUMABLES SELECTION GUIDELINES CHART

Abrasive Cutting Fluids		
Pace Product Name	Catalog Number	Packaging
MAXCUT [™] Abrasive Cutting Fluid (32 oz)	MAXCUT 1000-32	32 oz
MAXCUT [™] Abrasive Cutting Fluid (128 oz)	MAXCUT 1000-128	1 gallon
MAXCUT ~ 2 Cutting Fluid with Anticorrosion Additive (32 oz)	MAXCUT2 1000-32	32 oz
MAXCUT [™] 2 Cutting Fluid with Anticorrosion Additive (128 oz)	MAXCUT2 1000-128	1 gallon
300 mm (~12") MAXCUT [™] Abrasive Blades 3	32mm (~1¼") Arbor	
300 mm (~12") MAXCUT _™ Ahrasiye Blades 3	22mm (~1½") Arhor	
, ,	32mm (~1¼") Arbor	Packaging
Pace Product Name	` '	
Pace Product Name Soft Non-ferrous Materials (Aluminum, Brass, Zinc, Etc.)	Catalog Number	Packaging
300 mm (~12") MAXCUT™ Abrasive Blades 3 Pace Product Name Soft Non-ferrous Materials (Aluminum, Brass, Zinc, Etc.) Hard Non-ferrous Materials Soft Steels	Catalog Number MAX-E300	Packaging 10
Pace Product Name Soft Non-ferrous Materials (Aluminum, Brass, Zinc, Etc.) Hard Non-ferrous Materials	Catalog Number MAX-E300 MAX-C300	Packaging 10 10
Pace Product Name Soft Non-ferrous Materials (Aluminum, Brass, Zinc, Etc.) Hard Non-ferrous Materials Soft Steels	Catalog Number MAX-E300 MAX-C300 MAX-E300	Packaging 10 10 10



INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

4.5 Manual cutting (step-by-step procedure)

- 1. Connect electrical power and turn on master power switch on the side panel before opening the hood.
- 2. Position fixtures and sample near to cut-off wheel.
- 3. Close hood.
- 4. Turn on coolant and cutter motor.
- 5. Hand feed the cutting wheel or table feed into the sample.
- 6. Turn off cutter and allow blade to completely stop before opening cover. Note magnet safety switch has an approximately 5 second delay. Cutter hood cannot be opened until safety switch is deactivated.

4.6 Fixturing sample

For proper clamping, use the appropriate clamping vises to securely hold the sample in place. It is recommended that both sides of the part be clamped to avoid pinching of the blade (possibly breaking the blade) and to minimize burning of the workpiece during cutting.

For proper fixturing, take into account the initial stress on the samples.

Fixture examples



Fast clamping

Left lever (Cat. No. QCL-1100) Right lever (Cat. No. QCR-1100)

Note: vises have a slot in them so that the face does not rotate





UNIVERSAL Clamping vise

UNIVERSAL (Cat. No. MG15-01) Clamping height 0 - 3.5 in (0 - 90 mm)



INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

5.0 Maintenance

5.1 Introduction

The **MEGA-T300** requires very minimal maintenance. However, to increase the life of the saw, it is suggested that the cutting chamber be rinsed after cutting.

It is also recommended that the hood be left open in order to minimize corrosion inside of the chamber due to high humidity from the cutting operation and cutting fluid.

5.2 Cleaning outside cabinet

The cabinet and front shield should be cleaned occasionally with a moistened cloth. Do not use any chemicals or cleaning abrasives. We recommend an anti-condensation gel be used on the front panel.



INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

6.0 Trouble Shooting

More extensive trouble shooting, repair guides, video's, parts list are provided online at www.metallographic.com or

http://www.metallographic.com/PACE-service/MEGA-T300-service.html

Problem	Cause	Solution
No power or function	Unit is disconnected from main electrical power supply.	Verify electrical source and connection.
	b. Main power switch is off.	b. Turn on main power switch.
	c. Emergency stop button engaged.	c. Release by turning clockwise.
Main motor does not operate	Interlocking safety switch is in Unlock position (for shipping)	Turn safety switch to Lock position to operate cutter.
	b. Hood not closed.	b. Close hood
	c. Overload relays activated	c. Restart after resetting Q1 relay.
Pump motor does not operate	a. Plug is disconnected	a. Check plug.
Operate	b. Overload relay activated	b. Restart after resetting Q2 relay.
Excessive vibration during cutting	a. Not cutting with enough force	a. Increase cutting force.
during cutting	b. Specimen locked too far back on cutting table	b. Move specimen forward.
	c. Incorrect blade selection	c. Use correct blade.

Shipped in the Unlock position

Turn to Lock position to operate saw





If motor or pump does not turn on check Q1 - motor and Q2- pump relays

INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

----- **A**

7.0 List of Spare Parts

Part no.	Description	Images
	Electrical	
М3-СВ	MEGA-T300 Full control board replacement for MEGA T300	
GV2ME20C	QF1 overload switch (Schneider GV2-ME20C/13 -18A)	
GV2ME14C	MEGA-T300 QF1 overload switch (380V and 480V motor) (Schneider GV2-ME14C/6-10A)	Quintanii (
GV2ME05C	MEGA M250 QF2 overload switch (pump) (Schneider GV2-ME05C/0.63-1A)	The second secon
MEG-RCMU	RCMU timer (blade)	ON COLORES 3
MEG-RAMU	RAMU timer (safety switch)	3

Replacement of spare parts is possible with minimal expertise. Spare parts for recirculation system are also available upon request.

DACE TECHNOLOGIES

MEGA-T300 Abrasive Cutter

INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

Part no.	Description	Images
	Electrical	
LC1D18B7C	KM1, KM3 AC contactor (MEGA-T300) (Schneider LC1D12B7C)	
LC1D18B7C	KM2 AC contactor (MEGA-T300) (Schneider LC1D09B7C)	Pitchia
М3-F-Н	Fuse Holder (each)	Parison of the last of the las
M3-F1	MEGA-T300 F1 fuse (2A)	
M3-F2	MEGA-T300 F1 fuse (2A)	
M3-F3	MEGA-T300 F1 fuse (20A)	PACE
M3-F4	MEGA-T300 F1 fuse (20A)	
M3-F5	MEGA-T300 F1 fuse (2A)	
M3-F6	MEGA-T300 F1 fuse (2A)	
M3-MOTOR	MEGA-T300 Motor (4 KW)	PACE LOGIES TECHNOLOGICATION TECHNOLOGICATION TO THE PROJECT OF TH
M2-SAFE-T	MEGA M250 Safety lock out switch (top piece)	Quindra Control of the Control of th

INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

Part no.	Description	Images
	Electrical	
M2-SAFE-B	MEGA M250 Safety lock out switch (bottom piece)	PACE OCES TECHNOLOGICANO TECHNOLOGIC
M2-SAFE-K	Safety key	O SEPTIMENT
MEG-RL1	RL1 switch (24 V)	Control of the contro
RELAY-B	24V relay base	Bargreeter © C C
GTB	Ground terminal block	
ТВ	Terminal block	Pack Joe ES
MEG-LITE	Light housing	Participation of the second of

DACE TECHNOLOGIES

MEGA-T300 Abrasive Cutter

INSTRUCTION MANUAL

Part no.	Description	Images
	Electrical	
M-250B	Light blub (OSRAM Halogen DC24V 20 Watt)	C the control of the
M-SEAL	Wire plastic seal connector	And the second s
MEG-PUMP-220	Pump (220/380V)	Pack of the state
MEG-PUMP-480	Pump (400-480V)	
MEG-FPC	Female pump connector	PACE LOGIES TECHNIOS COMIC COM
MEG-MPC	Male pump connectors	PACE OCIES
MEG-AM1	Ammeter set-up kit	
MEG-R	Rectifier	

INSTRUCTION MANUAL

Part no.	Description	Images
	Electrical	
MEG-AM2	Amp overload light	
MEG-AM3	Amp overload light controller	A Bing Don't
MEG-TC-220	220-380V Transformer	ACE OGIETA
MEG-TC-480	440-480V Transformer	CE CLOPHICO
MEG-PS	Main power switch	

INSTRUCTION MANUAL

Part no.	Description	Images
	Electrical	
MEG-ON	On button (green)	CUTTER
MEG-LS	Light ON/Off switch	OFF ON PRODUCTION OF THE PRODU
MEG-LB	Power on light button	POWER
MEG-OFF	Off button (red)	Paraments STOP
MEG-W	Water ON/OFF Switch	AUTO MANUAL PUN
P150-210	Emergency Stop	And the same of th

INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

Part no.	Description	Images
	Mechanical-motor	
M3-SH	MEGA-T300 Blade shroud	Parista da de la companya del companya de la companya del companya de la companya
T3-M-001	Arbor bolt	
T3-M-002	Arbor inside flange	Questo de la constante de la c
T3-M-003	Arbor blade outside flange	
T3-M-004-1	Flange-lock-housing	Per de la contraction de la co
T3-M-004-2	Flange lock nut	
T3-M-004-3	Flange lock shaft	



INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

Part no.	Description	Images
	Mechanical-motor	
T3-M-003	Motor table base	Particular
T3-M-005	Motor table lock slot	Parit de come
T3-M-006	Motor base adjustment guide	Attach
T3-M-007	Motor base adjustment block	
T3-M-008	Motor extension shaft	A. C.
T3-M-009	Motor extension shaft housing	Acadacate



INSTRUCTION MANUAL

Part no.	Description	Images
	Mechanical-motor	
T3-M-010	Motor shaft	
T3-M-011	Motor shaft screw holder	
GS-T300	MEGA-T300 (400 Newton) gas cylinder pistons for hood	Andreas Andrea
T3-M-0012	Motor feed handle	Pace Cales
T3-M-004	MEGA-T300 blade locking rod	A Continue of the Continue of



INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

Part no.	Description	Images
	Mechanical-table	
M3-TG	MEGA-T300 Table feed guide rods (pair)	
М3-ВООТ	MEGA-T300 Table feed boot	Pace of other inches
M3-TFW	MEGA-T300 Table feed wheel	Pice to delegate the second delegate the secon
T3-M-013	MEGA-T300 Table feed screw	P. S.
T3-M-014	MEGA-T300 T-slot table	Pictorous .
T3-M-015	MEGA-T300 wheel feed locking lever	
T3-M-016	Table feed rail guide plate	. 25.

INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

Part no.	Description	Images
	Mechanical-table	
T3-M-017	Table feed back screw cover	
T3-M-018	Table feed rail guide block	
T3-M-019	Table feed sliding rail guide	
T3-M-020	Screw feed casting base outside cover-wheel	
T3-M-021	Screw feed casting base inside cover	
T3-M-054	Screw feed drive outside guide for table	
T3-M-022	Screw feed drive inside guide for table	Quita and a second
T3-M-023	Screw feed table wheel cover	C. C



INSTRUCTION MANUAL

Part no.	Description	Images
	Mechanical-table	
T3-M-024	Screw feed wheel-base cover-1	
T3-M-025	Table feed support bar large	
T3-M-026	Table feed support bar small	

Part no.	Description	Images
	Mechanical-tank	
MEG-TANK	Recirculating tank	PACE OLOGIES TECHNICION TOPHE COM TECHNICION TOPHE COM TECHNICION TOPHE COM TOPHE TOPHE TOPHE TOPHE COM TOPHE TOPHE TOPHE TOPHE TOPHE COM TOPHE
MEG-TANK-C1	Recirculating tank cover with square screen	PACE DIO GIES



INSTRUCTION MANUAL

Part no.	Description	Images
	Mechanical-table	
MEG-TANK-C2	Recirculating tank cover with round screen	PACE LOGIES TESTANDAMENTAL TOTAL CONTROL OF THE STATE OF
MEG-SCREEN- COV1	Tank screen/ single drain cover (MEGA-M250 / T300)	PACE LOGIES TECHNICATE AND THE COMMENT OF THE COMME
MEG-SCREEN-SQ	Coarse screen for MEGA abrasive cutter tanks (square)	Packologies Technologies
MEG-SCREEN-RD	Coarse screen for MEGA abrasive cutter tanks (round)	ACE OCE LE

DACE TECHNOLOGIES

MEGA-T300 Abrasive Cutter

INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

Part no.	Description	Images
	Mechanical-base	
M3-FEET	MEGA-T300 support feet (4/set)	J. Pring Dones
T3-M-028	MEGA-T300 lifting handles (2/set)	
T3-M-027	MEGA-T300 casting base	

Part no.	Description	Images
	Mechanical-Hood	
T3-M-029	MEGA-T300 Control panel fabricated cover	
T3-M-033	MEGA-T300 Cutting chamber hood	PACHOLOGIES C TECHNOLOGIES C



INSTRUCTION MANUAL

Part no.	Description	Images
	Mechanical-Hood	
T3-M-045	MEGA-T300 Right hood	Parting and the second
T3-M-030	Splash panel	ACE OLOGIES ECHNOLOGIES TECHNOLOGIES TOTAL CONTROL OF THE CONTR
T3-M-031	Water shield	
T3-M-032-M	Side panel (metal)	Picture of the second
T3-M-032-P	Side panel (plastic)	TECHNOLOGIES TECHNOLOGIES TECHNOLOGIES TECHNOLOGIES TECHNOLOGIES
МЗ-ВТ	MEGA-T300 rubber boot or side access	
М3-Т	MEGA-T300 template	Pace a cite Econolis



INSTRUCTION MANUAL

Part no.	Description	Images
	Mechanical-Hood	
T3-M-034	Cooling sheet top panel	Q. C.
T3-M-052	Electrical panel base plate	Q. in the second
T3-M-053	Electrical panel support box	Q
T3-M-035	Back sheet	
T3-M-036	Cooling back panel	Q. Harrier
T3-M-037	Side cooling vent	

INSTRUCTION MANUAL

Part no.	Description	Images
	Mechanical-Hood	
T3-M-038	Hood handle	
T3-M-039	Left front panel	Acceptages Rechifford Court of the Court of
T3-M-040	Lock shaft	
T3-M-041	Lock shaft block	
T3-M-043	Lock shaft ring	Q. Sandar
T3-M-044	M8 gas spring back holder	
T3-M-046	Right hood sheet cover	

INSTRUCTION MANUAL

Part no.	Description	Images	
	Mechanical-Hood		
T3-M-047	Rubber shield holder		
M3-FW	MEGA-T300 front safety window		
T3-M-048	Safety glass frame	Per de la contra del l	
T3-M-049	Small water shield		
T3-M-050	Water nozzle bracket	Charles and the second	
T3-M-051	Water split block		

INSTRUCTION MANUAL

Part no.	Description	Images
	Plumbing	
М3-D	MEGA-T300 Drain hose (40 mm) - 1.5m	DACE TECHNOLOGIES TECHNOLOGIES TECHNOLOGIES
М3-Н	MEGA-T300 Coolant line (20 mm) - 1.5 meter	PACE LOGIES TECHNOLOGIES THE CHARGE OF THE C
M3-COOL	MEGA-T300 Flexible coolant line (pair) (Part No. M3-COOL)	
MEG-RHV	Rinse hose and valve	P. British Control
MEG-RH	Rinse hose	Pack ordanes
M2-NOZ	Rinse hose nozzle	TECHNOLOGIES TECHNOLOGIES TECHNOLOGIES TECHNOLOGIES
M2-NOZ-T	Rinse hose nozzle tip	PACE COMES TECHNOLOGICS

INSTRUCTION MANUAL

Part no.	Description	Images
	Plumbing	
MEG-V	Rinse hose shut-off valve	PACE OGIES TECHNOLOGIANIC COM TECHNOLOGIAN PORTER METAL
MEG-IB	Brass inlet hose barb	
MEG-OB	Brass drain hose barb	

Part no.	Description	Images
	Tools	
G-GUN	Grease gun	PACE LOG
G-GUN-R	Grease refill	e e e e e e e e e e e e e e e e e e e
MEG-SQ	Square	PACE OLOGIES TECHNOLOGIES

MEGA-T300 Abrasive Cutter

INSTRUCTION MANUAL

Part no.	Description	Images
	Tools	
MEG-TB	Tool box	PACHOL GIES TECHNIQUE GOM
MEG-25HW	2.5 mm hex wrench	4 Patrickia positra,
MEG-10HW	M10 hex wrench	The standards
мз-wсн	MEGA T300 Wrenches (30/32 mm)	>= <

Part no.	Description	Images	
	Vice Components		
QC-REPAIR	Vice repair kit MEGA-T400 (kit include cam shaft and bar, spring and adjustment set screw for one vise)	Sign of the state	
MIG-01	Vertical clamping vise with 45 mm arm (1.8-inch) and shoe	Dac Moroculus m	
MG-15-01	Vertical clamping vise with 75 mm arm (3-inch) and shoe for MEGA-T400	Packarous de la constant de la const	




INSTRUCTION MANUAL

Part no.	Description	Images
	Vice Components	
МG-Н	65 mm (2.6-inch) height adapter for either MIG- 01 or MG15-01 vertical clamping vises	Action of the Comment
MG-TB	Vertical clamping vise t-slot bolt and nut	
MG-VW	Vertical clamping vise adjustment t-wrench	PACE OGIES TECHTOLOGICALINI
MEG-TS	T-slot nuts (2/set)	ace of the control of
QCR-1100	Right quick clamping vise for MEGA-T300 abrasive cutter	Pacificación de la companya de la co
QCL-1100	Left quick clamping vise for MEGA-T300 abrasive cutter	Pace a contest
MG-0159	Mechanical stop for cutting repetitive samples	PACE TECHNOLOGIES

INSTRUCTION MANUAL


3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

8.0 Electrical Drawing



MEGA-T300 Abrasive Cutter

INSTRUCTION MANUAL



MEGA-T300 Abrasive Cutter


INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

ITEM NO.	PART NUMBER	DESCRIPTION	QTY.
1	T3-M-033	Cutting chamber hood	1
2	GS-T300	Hood gas spring	2
3	T3-M-045	Right hood	1
4	T3-M-029	Control panel fabricated cover	1
5	T3-M-060-1	Champer hood hinge	1
6	T3-M-052	Electrical panel base plate	1
7	T3-M-053	Electrical panel support box	1
8	T3-M-037	Side cooling vent	1
9	T3-M-036	Cooling back panel	1
10	T3-M-034	Cooling top panel	1
11	T3-M-028	Liffting handle	4
12	T3-M-057	Three ways water splitter	1
13	M3-FEET	Support feet	4
14	T3-M-032	Splash panel	1
15	T3-M-038	Hood handle	1
16	M-D-45	45mm hose barb	1
17	M-D-20	20mm hose barb	3
18	T300-EX-A	Base-cover assembly	1

EX-Assembly

INSTRUCTION MANUAL






INSTRUCTION MANUAL

ITEM NO.	PART NUMBER	DESCRIPTION	QTY.
1	T3-M-027	Casting base	1
2	T3-M-011	Motor shaft screw holder	2
3	T3-M-058	Motor table holder	2
4	T3-M-010	Motor shaft	1
5	T3-M-003	Motor table base	1
6	T3-M-006	Motor base adjustment guide	1
7	T3-M-041	Lock shaft block	1
13	T3-M-023	Locking bracket	1
9	T3-M-040	Lock shaft	1
10	T3-M-043	Lock shaft ring	1
11	T3-M-015	Feed locking level	1
12	T3-M-062-1	Motor gas spring support	1
13	T3-M-005	Motor table lock slot	1
14	M3-MOTOR	Motor	1
15	T3-M-008	Motor extension shaft	1
16	T3-M-009	Motor extension shaft housing	1
17	M3-SH	Blade shroud	1
18	T3-M-0012	Motor feed handle	1
19	T3-M-002	Arbor inside flange	1
20	T3-M-049	Small water shield	1
21	T3-M-004-3	Blade locking rod	1
22	T3-M-030	Splash panel	1
23	T3-M-Bearing	Blade holder bearing	1
24	T3-10-12	12-inch blade	1
25	T3-M-051	Four ways water splitter	1
26	T3-M-042	Arbor blade outside flange	1
27	T300-IN-A	Base-Motor assembly	1

IN-Assembly

INSTRUCTION MANUAL



MEGA-T300 Abrasive Cutter

INSTRUCTION MANUAL

3601 E. 34th St. Tucson, AZ 85713 USA Tel. 520-882-6598 Fax 520-882-6599 email: pace@metallographic.com Web: http://www.metallographic.com

ITEM NO.	PART NUMBER	DESCRIPTION	QTY.
1	T3-M-014	T-slot table	1
2	T3-M-017	Table feed back screw cover	1
3	T3-M-021	Screw feed casting base inside cover	1
4	T3-M-054	Screw feed drive outside	1
5	T3-M-022	Screw feed drive inside	1
6	T3-M-013	Table feed screw	1
7	T3-M-020	Screw feed casting base outside cover	1
8	T3-M-025	Table feed support bar large	5
9	T3-M-026	Table feed support bar small	2
10	T3-QC-B	Quick claimping block	2
11	MG-TB	T-slot nuts	8
12	T3-M-019	Table feed sliding rail guide	2
13	T3-M-018	Table feed rail guide block	4
14	T3-M-050	Nozzle bracket	2
15	T3-M-024	Screw feed wheel-base cover	1
16	T3-M-016	Table feed rail guide plate	4
17	QCL-1100		1
18	QCR-1100		1

TB-Assembly